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COLTENE

Definition

DIATECH Z-Rex are diamond burs for high performance ceramics as for example zirconia and aluminum oxide.

Indications

- Endo access
- Crown adjustment
- Crown cutting

Recommendations for Use and Safe Operation in Dental Practices

The diamond instruments have to be selected (shape, size, type) according to the type of preparation. The ergonomic principles must be taken into account during the working procedures. The high-speed turbine or contra-angle handpiece as well as rotary instruments must be in perfect working condition. The instruments must be fully inserted and gripped in the high-speed turbine or contra-angle handpiece with utmost care. An optimal cutting performance and cutting life will be achieved by choosing the correct speed and pressure (depending on the substrate).

Before touching the preparation site with the instrument, make sure that the bur is rotating at the optimum speed (see packaging) and that the cooling system is working properly. Once the preparation has been completed, remove the rotary instrument from the site and allow it to come to a standstill.

The recommended force applied (reference value 1.5 N) is based upon the diameter of the instrument, the type of handpiece and the substrate to be worked on. The recommended speed is defined by the borderline circumferential velocities of 10 and 20 m/s. Use lower speed and pressure for finer grit size or for finishing, and a higher speed and pressure for coarse grit sizes and for bulk reduction. Jamming of the instrument can cause breakage. To prevent harm to the tooth structure, pulp and adjacent fillings, cool with sufficient amounts of water (minimum 50ml/min).

Whenever possible, the diamonds should be stored in their original packaging to enable them to be identified and traced. Instruments which are blunt, damaged, bent or no longer concentric, must be removed immediately and disposed of carefully. Only use flawless diamond instruments.

Recommended hygiene procedures in dental practices

It is recommended to sterilize diamond instruments prior to the first use. After use, instruments must be immersed immediately in a tank containing a suitable cleaning agent/disinfectant. Use a disinfectant/cleaning agent containing a corrosion inhibitor. The cleaning agent/disinfectant manufacturer's instructions must be obeyed. Strongly acidic and alkaline disinfectants may adversely affect the instruments.

If necessary, clean by mechanical means in advance, check for any remains and, if necessary, clean again. The instruments may be cleaned ultrasonically. In order to avoid the transmission of infectious germs, steam sterilization is recommended. Dry heat sterilization over 140°C damages the bur. Sterilizers must be used in accordance with the manufacturer's in-

structions.

Recommended speeds

Head diameter in 1/10 mm	Speed range (RPM)
012	160 000 – 300 000
014	135 000 – 275 000
016	120 000 – 240 000
018	105 000 – 210 000
023	85 000 – 165 000

Marking:

The  number is shown on the package.

Date of issue:

08-2017

Only supplied to dentists and dental laboratories or upon their instructions.

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COLTENE/WHALEDENT AG