Harvard BioCal® - Cem

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Harvard BioCal® - Cem

Bioactive, self-curing, self-adhesive resin reinforced MTA luting cement

Harvard BioCal® - Cem is a bioactive alkaline dual cure resin reinforced luting cement with MTA-fillers for permanent cementations.

Harvard BioCal® - Cem shows self adhesive properties, therefore application of an adhesive is not absolutely necessary. Nevertheless, application of a dual cure adhesive (e.g. Harvard Bond SE Dual) improves the bond strength of Harvard BioCal® - Cem to dentine and enamel and so guarantees a more reliable bond between tooth structure and Harvard BioCal® - Cem.

Indications

Permanent cementing of:

- Posts made of ceramic, metal and fibre-reinforced variants
- Crowns and bridges made of ceramic, zirconia and metal

Contraindications

Thin all-ceramic restorations of feldspathic ceramic or other low fusing ceramics, veneers, Maryland-bridges.

Side effects

In singular cases, the material may cause a sensitizing reaction in patients with a hypersensitivity to any of the ingredients. In these cases, the material should not be used.

Irritations resulting from direct contact with the pulp cannot be ruled out. Therefore for pulp protection areas close to the pulp should be covered with a thin layer of pulp capping material.

Incompatibility with other materials

Do not use in combination with substances containing eugenol because eugenol inhibits the polymerization of the material. Neither store the material in proximity of eugenol containing products, nor let the material come into contact with materials containing eugenol.

Preparing the MINIMIX-Syringe

Remove the cap of the MINIMIX-syringe and throw it away (do not use it again!). It is replaced by a supplied 1:1 mixing cannula. Turn the cannula 90° until it locks in position. The material is now ready for application.

Note:

Store used syringe with fixed used mixing cannula. Discard the first 2-3 mm of the extruded material. This has to be done for each new mix.

Store used syringe with fixed used mixing cannula in the dark.

The working time (23°C (74°F)) in the self cure mode is 1:30 minutes from start of mixing.

Application

Avoid water and saliva contamination during application and setting of the cement. For a dry working area adequate isolation of the tooth is required.

1. Cementation of posts

1.1. Preparing of the root canal

Use of a rubber dam to isolate the tooth is strongly recommended.

Prepare the selected post according to manufacturer directions.

Treat the root canal endodontically and clean it in the usual manner (e.g. with sodium hypochlorite solution, approx. 3%). Rinse with water and dry with paper points.

Optional:

Apply a suitable dual cure bonding agent (e.g. Harvard Bond SE Dual) according to the corresponding instructions.

1.2. Filling of the root canal

Prepare the MINIMIX-syringe as described above and attach an Endo-Tip to the mixing cannula. Insert the Endo-Tip as deep as possible into the root canal and apply Harvard BioCal® - Cem starting apical. Keep the end of the Endo-Tip in the material to avoid bubbles. When the root canal has been completely filled remove the Endo-Tip from the material.

If necessary apply a thin layer of Harvard BioCal® - Cem on the post. Place the post in the filled root canal and hold in position under moderate pressure. Light cure the coronal part of the cemented post for **20 seconds** with a polymerization unit (wavelength range of 400-500 nm and light intensity of at least 1000mW/cm²).

2. Cementation of crowns and bridges

2.1. Preparing

Prepare the luting side areas of the restoration with a suitable bonding agent according to the corresponding instructions:

2.2. Preparation of cavity/tooth stump

For a sufficient retention the height of the prepared stump must be at least 4 mm with a steep preparation angle (maximum 6 degrees).

Prior to final cementation, clean the prepared stump or cavity thoroughly with pumice slurry, rinse with water and slightly dry in 2-3 intervals with oil free air. **Do not overdry!**

Optional:

Apply a suitable dual cure bonding agent (e.g. Harvard Bond SE Dual) according to the corresponding instructions.

2.3. Cementing of crowns and bridges

For cementing apply a uniform coat (0.5 mm) of Harvard BioCal® - Cem on the luting side of the crowns and bridges and on the prepared tooth areas. Seat the restoration under light pressure and retain it in place until the cement has completely cured. Remove excess material after brief light exposure (approx. 5 sec).

Light cure all marginal areas of the restoration from each direction (mesio-oral, disto-oral, mesio-buccal, disto-buccal) for **20 seconds** with a polymerization unit (wavelength range of 400-500 nm) and light intensity of at least 1000mW/cm².

The dual cure cement system will auto cure within **4 minutes**.

Additional notes

- Unpolymerized material may have an irritating effect and may lead to a sensitizing reaction against methacrylates.
- Avoid contact with skin, mucous membrane and eyes.
- If the material comes into contact with skin, immediately wash with water and soap.
- Contact with eyes may cause severe eye damage. Wear eye protection. If the material comes into contact with eyes, immediately rinse with copious amounts of water and seek medical advice if required.
- Commercial medical gloves do not protect against the sensitizing effect of methacrylates.

Storage

Store at 2 – 8 °C (36 – 46 °F) in the original seal-pack. After opening of the seal-pack store the material protected from light and moisture at room temperature (19 – 25 °C / 66 - 77 °F). Use up within 6 months and before the end of expiry date. Do not use after expiry date.

Warranty

Harvard Dental International GmbH warrants this product will be free from defects in material and manufacture. Harvard Dental International GmbH makes no other warranties including any implied warranty of merchantability or fitness for a particular purpose. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusive remedy and Harvard Dental International GmbH's sole obligation shall be repair or replacement of the Harvard Dental International GmbH product.

Limitation of liability

Except where prohibited by law, Harvard Dental International GmbH will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

For dental use only!

Store product out of reach of children!

Keep the instructions for use for the duration of the application.

Order-No.

7100011

Article

Harvard Harvard BioCal® - Cem
2.5 ml minimix syringe, 10 mixing tips

 Manufacturer	Harvard Dental International GmbH Margaretenstr. 2 – 4 15366 Hoppegarten, Germany	 0482 Made in Germany	Tel: +49 (0) 30/99 28 978-0 Fax: +49 (0) 30/99 28 978-19 info@harvard-dental-international.de www.harvard-dental.de	V03 2022-03
---	--	--	--	-------------